

# Dragon Medical One Quick Tips

## Useful commands

- New line/new paragraph
- Delete that/Scratch that (deletes the last dictated word(s))
- Select <word> and then delete that (to delete a specific word or phrase)
- Select <word>
- Select <word> through <word>
- Select that (selects last utterance)
- Select all
- Insert [before/after] <word>
- Undo that
- Redo that

## Navigating

- Go to [top/bottom/end]
- Go to [start/end] of [field/sentence/paragraph]
- Go back (returns the cursor to the previous location after using a command)
- Resume with <word> (deletes all text between <word> and the current cursor position)

## Getting help

- Say "Open help" or "What can I say?"
- Say "Open training" to launch short instructional Dragon Medical One videos

## Adding words to your vocabulary from dictated text

- Say "Select <word(s) to add>"
- Say "Add word." Click Use Default Pronunciation or click the microphone icon and say the word/phrase

## Microphone management\*

Mic On .....		
Mic Off .....		
Mic in Standby Mode .....		

\*Your DragonBar may look different depending on your version of Dragon Medical One.

**For more quick tips, just say "Open training" from the Dragon Medical One application.**

MEDICAL DEVICE DISCLAIMER. Nuance products are not designed, intended or made available as a medical device(s) and are not designed or intended to be a substitute for professional medical advice, diagnosis, treatment, or judgment and should not be used to replace or as a substitute for professional medical advice, diagnosis, treatment, or judgment.

Do Not Distribute. Contains Nuance Information shared under NDA.