

# High-quality dictation microphone with mouse functions.

Boost clinician productivity with Nuance PowerMic 4 ergonomic control of dictation functions and on-screen field navigation.

Nuance PowerMic 4\* brings clinician dictation to a new level by making it quick and easy to control the microphone, navigate and select on-screen fields, and insert, review, and edit generated text.

## Control at your fingertips

With new features such as a larger pointing device, additional buttons, and magnetic mount support, clinicians will enhance productivity with an improved ergonomic control of dictation and speech recognition functions.

With Nuance PowerMic 4, clinicians can rapidly navigate EHR, template, or report fields and insert text using dictation for new levels of usability and efficiency.

---

## KEY FEATURES

- **USB plug-and-play installation** ensures compatibility with a wide range of computers.
- **Unidirectional noise-cancellation microphone** ensures performance accuracy in even the noisiest environments.
- **Easy-to-clean surface** promotes safety in clinical settings.

---

## AVAILABLE MODEL

- 2,7 m USB cable
- 



\* Version requirements  
– Dragon Medical One version 2021.4 SR1  
– Dragon Medical SpeechKit version 2021.4

## Technical specifications

---

**Operating System:** Windows® 10, Windows® 11

---

**USB:** USB 2.0, type A (compatible with USB 3.0)

---

**Supported Environments:** Citrix® and VMware

For more information, please contact Nuance Technical Support.

---

### Frequency Response Input/Output

Microphone 20 - 12,000 Hz

Speaker 300 - 7,500 Hz

Signal-to-Noise Ratio: >70dbA

---

### Power Consumption

Current: < 100mA (max playback volume)

Power: max 0.5 VA (max playback volume)

---

### Operation Conditions

Operating Temperature 5° to 45° C (41°F to 113°)

Humidity: 10% - 90%

---

### Emissions – Class B

— EN 55032:2015 – Electromagnetic Compatibility of multimedia equipment—Emissions requirements

---

### Immunity

— EN 55035:2018 – Electromagnetic Compatibility of multimedia equipment—Immunity requirements

---

### Compliance

#### Manufacturers declaration of conformity

The device complies with the following EU Directives

- 93/68/EEC—CE marking directive
- 2014/30/EU—Electromagnetic Compatibility
- 2011/65/EU & 2015/863/EU—RoHS

Date of first affixing the CE mark 2021



## LEARN MORE

To find out more about the Nuance PowerMic 4 microphone please visit [nuance.com/powermic-no](https://nuance.com/powermic-no)



### About Nuance Communications, Inc.

[Nuance Communications](https://nuance.com) is a technology pioneer with market leadership in conversational AI and ambient intelligence. A full-service partner trusted by 77 percent of U.S. hospitals and 85 percent of the Fortune 100 companies worldwide, Nuance creates intuitive solutions that amplify people's ability to help others. Nuance is a Microsoft company.